
NeuWave™

Ablation Probes

en	Instructions for Use
fr	Mode d'emploi
de	Gebrauchsanleitung
es	Instrucciones de uso
it	Istruzioni per l'uso
nl	Gebruiksaanwijzing
sv	Bruksanvisning
da	Brugsvejledning
pt	Instruções de utilização
fi	Käyttöohjeet
ko	사용 지침
zh-t	使用說明

REF PR15

REF PR20

ETHICON

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

This Page Intentionally Left Blank

NEUWAVE™ Ablation Probes

Ablation Probes Instructions for Use

REF

PR15


NEUWAVE PR Probe
15CM 17ga

REF

PR20

NEUWAVE PR Probe
20CM 17ga

Read this entire document and any accessory documents prior to use.
Refer to the NEUWAVE™ Microwave Ablation System User Manual for complete instructions for use.

For use only with the NEUWAVE Microwave Ablation System.

NEUWAVE Ablation Probes are provided sterile. Follow facility procedures for sterile device handling.

CAUTION U.S. Federal law restricts this device to sale by or on the order of a physician. Outside the U.S.A., check local laws for any restrictions that may apply.


NEUWAVE is a registered trademark of NeuWave Medical, Inc.

Other brand names or product names used in this manual are trademarks or registered trademarks of their respective holders.

Symbols used in instructions for use or on probes

Symbols replace words on the probes, on the display, or in NeuWave manuals and Instructions for Use. The following symbols appear on the NEUWAVE Microwave Ablation Probes.

	Manufacturer
	Sterilized using Ethylene Oxide Gas
	The expiration date
	A single-use product that should not be reused
	Authorized representative in the European Community
	Type BF Equipment protection against electric shock
	Refer to user manual
	Device complies with the European Medical Device Directive (MDD)
	Device contains no natural rubber
	Caution: Federal law in the US restricts this device to sale by or on the order of a physician.
	Keep dry
	Do not use if package is damaged
	Storage humidity range
	Storage temperature range
	The lot number of the item, unique to a specific number of units manufactured under identical conditions.
	Is the reference or part number of a given item. REFs are identical for all items of a given type.
	Unit quantity

Indications For Use:

The NEUWAVE™ Microwave Ablation System and Accessories are indicated for the ablation (coagulation) of soft tissue in percutaneous, open surgical and in conjunction with laparoscopic surgical settings, including the partial or complete ablation of nonresectable liver tumors.

The NEUWAVE™ Microwave Ablation System is not indicated for use in cardiac procedures.

Probe disposal

Disconnect all ablation probes from the PDM (Power Distribution Module) and dispose probes per your facility's sharps and biohazard handling procedure.

Probe Use

- WARNINGS**
- The NEUWAVE Ablation System should not be used in the following situations:
 - Pregnant patients – potential risks to patient and/or fetus have not been established.
 - Patients with implantable pacemakers or other electronic implants. Implanted electronic devices may be adversely affected by microwave power.
 - Use on the central nervous system.
 - Electrosurgical/electrocautery devices may interfere with the NEUWAVE Microwave Ablation System and cause system errors. Ensure that all NEUWAVE probes are removed from the patient prior to using electrosurgical/electrocautery devices.
 - Use only NEUWAVE probes with the NEUWAVE Microwave Ablation System. Probes from other manufacturers may cause patient injury or fail to function properly.
 - This product should only be used by physicians and staff properly trained in the use of the technology and its associated warnings and cautions.
 - NEUWAVE Ablation Probes are provided sterile and are single use devices. Do not reuse or re-sterilize any product labeled “SINGLE USE”. Doing so may result in cross-contamination, injury to the patient or medical staff, or equipment malfunction.
 - Do not use damaged probes, probes whose sterile barrier is suspect and probes beyond the expiration date noted on the packaging.
 - Use medical imaging devices to verify proper probe placement prior to delivering microwave energy to the patient.
 - Intraprocedure and/or post-ablation imaging are recommended to assess the extent of tissue coagulation.
 - Probe tips are sharp. Handle with care.
 - Always inspect probes and cables for chips, cracks, or other damage before each use. Do not use damaged probes.
 - Do not defibrillate a patient with a probe inserted. Completely remove the probe from the patient before defibrillation.
 - Prior to moving patient into a CT imager, ensure the placed probes have sufficient clearance from the CT bore. Patient injury can result if the probe handles contact the CT bore.
 - Do not place probes so that the probe cables are draped on or near the patient's head or neck.
 - Prior to starting an ablation, use imaging to confirm proper probe placement and that the probe is not bent or broken.
 - Do not attempt to bend or reshape probes. Probe malfunction may occur, possibly causing patient injury.
 - Probes are provided sterilized. Follow your facility sterile handling guidelines.
 - The probe cable can heat and reach 60° C during energy delivery. When placing the probe, verify that the probe cable does not rest on the patient's skin. Use the included clips to secure the probe cable away from the patient's skin as needed.
 - The non-active region of the probe cannula can heat and reach 48° C during energy delivery.
 - Any undue handling or touching of the probe shaft during or following use could result in a thermal injury to the patient or user.
 - Never bend or apply excessive force to the probe. Probe malfunction may occur possibly causing patient injury.
 - Never remove a probe from the patient while “Tissu-loc” mode is enabled or the displayed temperature is below 0° C. Serious internal injuries may result.
 - Never activate a probe function other than “Test” when the probe is not placed in tissue.
 - Probes removed from the body may be hot. Handle with care to avoid user and/or patient injury.
 - Do not connect ablation probes directly to the back of the cart as potential injury to the patient or user is possible.
 - Do not attempt to disconnect a probe when any functions are active.
 - Do not handle probes that appear excessively cold.
 - NEUWAVE Ablation Probes can be used with needle introducers. If a needle introducer is used, it must be retracted so that it doesn't interfere with the planned ablation zone prior to delivering energy. Failure to retract the needle introducer a sufficient distance may impact energy delivery.
 - Do not apply a clamp or hemostat to the probe shaft or cable as this may damage the probe and cause a malfunction.
 - Probe movement during ablation is possible. Probe movement can be caused by patient movements from breathing, coughing, etc. and also by pressure applied to the probe when target tissue contracts in response to the ablation. To help prevent probe movement, hold the probe handle in place during the initial period of the ablation (approximately 45 seconds). Monitor the probe for movement throughout the procedure and hold the probe in place as needed. If repositioning of a probe is required, stop energy delivery prior to repositioning the probe.
 - Charred tissue may accumulate on the probe tip during planar coagulation. As needed, use a soft sterile soft cloth or soft pad to remove any tissue build up on the probe tip. Probe tip may be hot. Never use a rough surface or scratch pad to clean the probe tip as this may damage the probe and result in errors or malfunction.
 - There are temperature sensors on the outside of the probe shaft. If the temperature sensors are damaged during placement or use of the probe, the NEUWAVE Microwave Ablation System will generate an error and disable the probe. To minimize the risk of damaging the temperature sensors during placement, avoid penetrating rigid structures such as bones and cartilage without the use of an introducer. Sharp instruments and hemostats/clamps should never be used along the probe shaft. Ultrasound guides and needle introducers should be used with caution.
 - The Surgical Clip (DR-001193/RCPK) is the only clip that is designed for use with 17-gauge and 15-gauge probes. To avoid damage to probes, do not use probes larger than 17-gauge with the Surgical Clip (DR-000387/RC17PK).
 - Keep metal objects, including other medical instruments away from the probe shaft when delivering energy. Energy from the probe may radiate to the metal objects causing them to heat.
 - Ensure the tip of the probe is inserted in tissue to the minimum insertion depth marking prior to activating power.
 - When performing procedures, always use the lowest power setting and shortest time that will achieve the intended results. Refer to Ex-Vivo Ablation Results for example ex-vivo ablation sizes based on power and time settings.

CAUTION

Important

- Probe temperatures are highly dependent on tissue variables such as perfusion and vasculature. Probe temperatures are likely to vary from procedure to procedure based on these tissue variables. These temperature variations do not, by themselves, indicate a malfunctioning probe or system.

Target Ablation Use Guide

Single Probe			
Lung Tissue	Liver Tissue	Kidney Tissue	Other Soft Tissue
✓	✓	✓	See 1 below
Multiple Probes			
Lung Tissue	Liver Tissue	Kidney Tissue	Other Soft Tissue
See 3 below	✓	See 1 and 2 below	See 1 below

- Performance characteristics of NEUWAVE Ablation Probes have not been established in all probe and tissue type combinations.
- Multiple probes make larger ablations than single probes. When using multiple probes, ensure that the target tissue is an appropriate size to limit the risk of unintended thermal damage to non-target tissue.
- Performance characteristics of multiple probes and 15ga or larger probes have not been established for lung tissue.

NEUWAVE PR probes were designed to limit the length of the ablation for instances when a shorter ablation zone is desired. NEUWAVE PR probes were developed to provide physicians with an additional ablation probe designed specifically for ablating smaller targets. The NEUWAVE PR probes are designed to produce ablations that quickly encompass the tip of the probe while limiting the overall length of the ablation.

In percutaneous applications, always use imaging to confirm desired ablation results. See tables of ex-vivo ablation results.

Up to three (3) probes can be connected to the NEUWAVE Microwave Ablation System and powered at any time. This can be done to create one, larger ablation or three (3) independent ablations. Refer to tables of ex-vivo ablation results and guidelines for minimum and maximum probe separation distances.

Refer to the NEUWAVE Ablation System User Manual for complete instructions for use.

A protective sheath is provided installed on the NEUWAVE Ablation Probes to protect the user from the sharp probe tip during removal from the package and to prevent accidental contact of the probe tip from non-sterilized surfaces during setup. The probe protective sheath must be removed prior to probe testing.

When using introducers during bone ablations, reference guidance in the NEUWAVE Microwave Ablation System User Manual.

Surgical Coagulation Use Guide

NEUWAVE 17-gauge probes can be used for Surgical Coagulation in either a single or dual probe configuration. When using two 17-gauge probes, NeuWave recommends using the Surgical Clip (DR-001193/RCPK). See the NEUWAVE Microwave Ablation System User Manual for more information.


Ex-Vivo Ablation Results

The following are examples of ex-vivo ablation results when using single and multiple probes at various power and time settings. These examples are from ex-vivo animal (bovine and porcine) tissue.

The size of ablation zones may vary in clinical settings due to vasculature structure and the cooling effect of blood flow not present in ex-vivo tissue.

Single Probe Ablations

Model NEUWAVE PR Probes (PR15 and PR20)

Bovine Liver – Ex-Vivo

Power	Time	Length (B)	Diameter (A)	Distance Past Tip (C)	Tip to Center (D)	Volume
30W	5 Minutes	2.6cm	1.8cm	0.3cm	1.0cm	4.7cm³
30W	10 Minutes	3.2cm	2.4cm	0.4cm	1.2cm	9.6cm³
65W	5 Minutes	4.0cm	2.5cm	0.5cm	1.5cm	13.6cm³
65W	10 Minutes	4.6cm	3.1cm	0.6cm	1.7cm	23.0cm³

Porcine Kidney – Ex-Vivo

Power	Time	Length (B)	Diameter (A)	Distance Past Tip (C)	Tip to Center (D)	Volume
30W	5 Minutes	2.6cm	1.8cm	0.4cm	0.9cm	4.3cm³
30W	10 Minutes	3.2cm	2.3cm	0.5cm	1.1cm	8.6cm³
65W	5 Minutes	3.9cm	2.6cm	0.5cm	1.4cm	13.5cm³
65W	10 Minutes	4.5cm	3.1cm	0.6cm	1.6cm	22.0cm³

Porcine Lung – Ex-Vivo

Power	Time	Length (B)	Diameter (A)	Distance Past Tip (C)	Tip to Center (D)	Volume
30W	5 Minutes	2.2cm	1.3cm	0.2cm	0.9cm	2.0cm³
30W	10 Minutes	2.7cm	1.6cm	0.3cm	1.1cm	3.8cm³
65W	5 Minutes	3.6cm	1.9cm	0.4cm	1.4cm	7.1cm³
65W	10 Minutes	4.1cm	2.3cm	0.4cm	1.6cm	11.0cm³

Multiple Probe Ablations

Based upon ex-vivo results, NeuWave recommends a maximum of 2 cm spacing between probes. This results in a contiguous ablation. Ablations spaced more than 2 cm had non-contiguous results. Ablations spaced a minimum of 8 cm resulted in 3 independent ablations zones that did not overlap.


Performance characteristics of multiple probes has not been established for lung tissue.

Multiple probes make larger ablations than single probes. When using multiple probes, ensure that the target tissue is an appropriate size to limit the risk of unintended thermal damage.

Model NEUWAVE PR Probes (PR15 and PR20)


2 Probe: 65 W Power Setting in Bovine Liver, Ex-Vivo

Time	Length (B)	Diameter (A)	Probe Spacing (S)	Volume
10 Minutes	5.8cm	4.6cm	2.0cm	63.1cm³

3 Probe: 65 W Power Setting in Bovine Liver, Ex-Vivo

Time	Length (B)	Diameter (A)	Probe Spacing (S)	Volume
10 Minutes	6.5cm	5.3cm	2.0cm	93.5cm³

Planar Coagulation Technique

The following results were obtained by following the technique for Planar Coagulation identified in the NEUWAVE Microwave Ablation System User Manual (for revision 3.1.X and higher) using a NEUWAVE PR Probe, a 0.75 cm probe spacing and coagulating in ex-vivo bovine liver and porcine kidney.

Model NEUWAVE PR Probes (PR15 and PR20)

Single Probe: In Bovine Liver – Ex-Vivo – Thermal Spread

Power	Time	Avg Depth	Avg Lateral
65W	20 sec	1.7cm	1.3cm
95W	20 sec	1.9cm	1.5cm

Dual Probe: In Bovine Liver – Ex-Vivo - Thermal Spread

Power	Time	Avg Depth	Avg Lateral
65W	20 sec	1.6cm	1.3cm
95W	20 sec	1.8cm	1.6cm

Single Probe: In Porcine Kidney – Ex-Vivo - Thermal Spread

Power	Time	Avg Depth	Avg Lateral
65W	20 sec	1.1cm	0.8cm
95W	20 sec	1.3cm	0.9cm

Ablation Zone Guidance Factors

The ablation zone data provided on-screen and in the instructions for use are from ablations performed in ex-vivo animal tissue and do not account for the cooling effects of blood flow.

Blood flow will generally result in a smaller ablation.

Lower power ablations are more impacted by the cooling effects of blood flow.

Additional factors that may impact the final ablation zone include but are not limited to:

Factors that may result in a smaller ablation	Factors that may result in a larger ablation
Increased/high level of vascularity	Decreased/low level of vascularity
Proximity to large blood vessels	Prior treatments such as TACE or embolization

Physicians should consider all factors when determining the appropriate ablation parameters for each case.

For Canadian Customers Only:

Contraindications:

The NEUWAVE™ Microwave Ablation System is contraindicated for use in the following situations:

- Pregnant patients – potential risks to patient and/or fetus have not been established.
- Patients with implantable pacemakers, neurostimulators, or other electronic implants. Implanted electronic devices may be adversely affected by microwave power.
- Use on the central nervous system.

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©2021 NeuWave Medical, Inc. All Rights Reserved.
Johnson & Johnson Medical GmbH